



Fluorochrome UniMF Microsphere

Description

Microsphere Composition: Melamine-Formaldehyde Resin

Fluorescent Composition:

1. Nile Blue A (red fluorescence; Excitation/Emission = 636/686nm)
2. Rhodamine B (orange fluorescence; Excitation/Emission = 550/584nm)
3. Rhodamine 6G

Form : Aqueous dispersion or Dry powder

Approximate Concentration (W/V): 5% or 10% solids

Sodium Azide Concentration : 50 ppm

Surfactant(W/V) : <0.1% or None

Physical Data

Nominal Diameter(μm) : 6, 8, 13, 50

Density (g/cm^3) : 1.51


Refractive Index: 1.68 (589nm, 25°C)

Physical and Chemical Properties


Hydrophilic surface, possess many functional groups on the surface (methylol groups, amino groups), which can be used for a covalent attachment of other ligands; High cross linking density, High temperature stability up to 300°C; Superior mechanical strength; Stable in acids and bases; Extremely high stability in organic solvents, no swelling or shrinking upon contact with organic solvents; Soluble in organic solvents like benzene, halogenated hydrocarbons, or acetone; Reduced non-specific protein binding activity.

Fluorescence Microscopy Image

Rhodamine B labeled MF microsphere




Green light excitation




Blue light excitation

Nile Blue A labeled MF microsphere




BRI (white light)




Green light excitation

Fluorescence Emission Spectra of Microsphere



Rhodamine B



Nile Blue A

Storage and Handling

Aqueous dispersions of melamine resin particles have excellent stability. Storage at room temperature is possible without bacterial growth. Particles can be washed with alcohol, air dried and autoclaved. Dried particles can be redispersed in water without any agglomeration. Dispersions of MF particles can be frozen.

Application

1. MF particles find wide applications as model systems in medicine, biochemistry, colloid chemistry, and aerosol research;

2. The fluorescent particles can be used as standards (e.g. in flow cytometry, confocal laser scanning microscopy, light scattering instruments) as well as tracers in environmental science, flow measurements in gases and liquids like Laser Doppler Anemometry (LDA), Particle Dynamics Analysis (PDA), Particle Image Velocimetry (PIV), Digital Imaging Velocimetry (DIV) and Laser Speckle Velocimetry (LSV); In applications with a high background light level, fluorescent seeding particles can significantly improve the quality of vector maps from PIV and LDA measurements.

2. 1 Single-phase liquid flows —The enhancement of image contrast is beneficial for most PIV applications even for simple liquid flows in order to improve the measurement accuracy.

2. 2 Multi-phase flows —Bubbly flows and droplet flows etc.

2. 3 Industrial large-scale flows —Mixing tank and filtering flows etc.

2. 4 Stereo PIV —The emission of fluorescence is independent on the viewing angle unlike the Mie scattering. The intensity difference often seen between cameras can be minimized, leading to enhance the accuracy of Stereo PIV.

2. 5 Micro PIV

2. 6 Near-wall flows